



# HIS Integration

## Product Specifications

### HL7

Supports HL7 version 2.x

**Inbound** ADT, SIU or ORM

**Outbound** Standard ORM

Custom inbound or outbound messages  
(Available from Olympus Professional Services)

### Network

Fully integrated with VaultStream Patient

**Network** Configurable (TCP/IP)

**Port**

### DICOM

**Inbound** Modality Work List (MWL)

Modality Performed Procedure Step

**Outbound** (MPPS), Visible Light (VL) image export;  
endoscopic or photographic (Video)

### Workflow

Automatic or pick-list patient registration and discharge

Support for update and cancellation of orders

Results Integration

Hyperlink to VaultStream

**Results access options** Link to PDF report file

HL7 encapsulated result PDF

Full user authentication to the VaultStream application  
through Active Directory

## Features


Integrate nCare with hospital applications for seamless workflow

Reduce manual interventions and errors

Streamline patient registrations

Provide a single point of access for images and patient data

Support unified patient record



Improved Clinician Efficiency with Orders and Results Integration

— Inbound/Order — Outbound/Results

— VNA — VNA — VNA — VNA